

# Next-generation test, improved report

Page 1: Quick, clear access to the information you need most<sup>1</sup>

**ABOUT THE TEST** FoundationOne® Liquid CDx is a next generation sequencing (NGS) assay that identifies clinically relevant genomic alterations in circulating cell-free DNA.

**PATIENT**  
 DISEASE Lung adenocarcinoma  
 NAME 999999999, FR  
 DATE OF BIRTH Not Given  
 SEX Not Given  
 MEDICAL RECORD # Not Given

**PHYSICIAN**  
 ORDERING PHYSICIAN Not Given  
 MEDICAL FACILITY Not Given  
 ADDITIONAL FACILITY Not Given  
 MEDICAL FACILITY ID Not Given  
 PATHOLOGIST Not Given

**SPECIMEN**  
 SPECIMEN ID Not Given  
 SPECIMEN TYPE Blood  
 DATE OF COLLECTION Not Given  
 SPECIMEN RECEIVED Not Given

PATIENT 999999999, FR TUMOR TYPE Lung adenocarcinoma REPORT DATE 28 August 2020  
 COUNTRY CODE FR ORDERED TEST # ORD-XXXXXX-XX

**Genomic Signatures**  
 Blood Tumor Mutational Burden - 5 Muts/Mb  
 Microsatellite status - Cannot Be Determined  
 Tumor Fraction - 13%

**Gene Alterations**  
 For a complete list of the genes assayed, please refer to the Appendix.  
 EGFR exon 19 deletion (L747\_A750>P)  
 TP53 R267P

5 Therapies Approved in the EU 10 Clinical Trials  
 0 Therapies with Lack of Response

**GENOMIC SIGNATURES**

**Blood Tumor Mutational Burden - 5 Muts/Mb**

**Microsatellite status - Cannot Be Determined**

**Tumor Fraction - 13%**

GENE ALTERATIONS	VAF %
EGFR - exon 19 deletion (L747_A750>P)	0.20%

10 Trials see p.11

**THERAPY AND CLINICAL TRIAL IMPLICATIONS**

No therapies or clinical trials. see Genomic Signatures section

Unable to determine Microsatellite status due to insufficient evidence of genomic instability.

Tumor fraction is an estimate of the percentage of circulating-tumor DNA (ctDNA) present in a cell-free DNA (cfDNA) sample based on observed aneuploid instability.

THERAPIES APPROVED IN THE EU (IN PATIENT'S TUMOR TYPE)	THERAPIES APPROVED IN THE EU (IN OTHER TUMOR TYPE)
Afatinib 1	
Dacomitinib 1	
Erlotinib 1	
Gefitinib 1	
Osimertinib 1	

NCCN Category

Electronically Signed by Julia A. Elvin, M.D., Ph.D. • 01 June 2020  
 Julia Elvin, M.D., Ph.D., Laboratory Director CLIA: 22D2027531  
 Shakti Ramkissoon, M.D., Ph.D., M. M. Sc, Laboratory Director CLIA: 34D2044309  
 Foundation Medicine, Inc. • 1-888-988-3639

Sample Preparation: 150 Second St., 1st Floor, Cambridge, MA 02141 • CLIA: 22D2027531  
 Sample Analysis: 150 Second St., 1st Floor, Cambridge, MA 02141 • CLIA: 22D2027531  
 Post-Sequencing Analysis: 150 Second St., 1st Floor, Cambridge, MA 02141 • CLIA: 22D2027531

PAGE 1 of 17

6

7

- 1 **Genomic signatures.** bTMB and MSI\* status, that may help inform eligibility for immunotherapy
- 2 **Gene alterations.** Clinically-relevant alterations in >300 tested cancer-related genes
- 3 **Pertinent negative results.** Rules out important alterations that are not present
- 4 **Therapies with clinical benefit.** Therapies approved in the EU\*\* for your patient's genomic signatures and gene alterations
- 5 **Clinical trials.** Relevant trials that your patient may be eligible for, based on their genomic profile and geographic location
- 6 **Gene alterations with no reportable options.** To help you rule out uncertainty and determine the most appropriate course of action
- 7 **Variant allele frequency percentage (VAF%).** for base substitutions and insertions

Reports vary according to regional differences, e.g. EU reports list EU-approved therapy options to support clinical decision-making.  
<sup>1</sup>FoundationOne Liquid CDx reports MSI-H status.  
<sup>\*\*</sup>Therapies contained in this report may have been approved through a centralised EU procedure or a national procedure in an EU Member State. For additional information on the NCCN categories please refer to the NCCN Compendium\* (www.nccn.org).  
 MSI, microsatellite instability, bTMB, blood tumour mutational burden.

# Page 2 onwards: detailed supporting information and guidance

Gives important background information on **genomic signatures and gene alterations** in your patient's tumour, notably the evidence behind associated potential treatment strategies.

Provides evidence-based insights on the therapies that **match your patient's profile, in their tumour type and others.**

Details **ongoing clinical trials** that your patient may be eligible for, with trial NCT number, phase, title and locations.

Lists **variants of unknown significance (VUS)**, included in the event that they become clinically meaningful in the future.

Tells you **what we tested** to give you a comprehensive picture of our service.

# FoundationOne Liquid CDx report: For insights for your patients at optimal times beneficial to their treatment journey

FOR MORE INFORMATION PLEASE CONTACT YOUR ROCHE REPRESENTATIVE

www.foundationmedicine.se | 08-726 11 00

1. Data on file: FoundationOne Liquid CDx mock report, 2020. Foundation Medicine® and FoundationOne® are registered trademarks of Foundation Medicine®, Inc. Roche is the licensed distributor of Foundation Medicine® products outside of the United States www.foundationmedicine.se | 08-726 11 00 M-SE-00000431 | Date of preparation: September 2020

